



germanengineering

BG 21 Glas-Diagnosewaage

- Grafisches Druckmotiv mit "Durchblick"
- Großes LCD-Display mit Multi-Anzeige
- Ziffernhöhe (Gewichtsanzeige): 25 mm
- Körpergewicht, -fett, -wasser, Muskelanteil, Knochenmasse, Kalorienbedarf (BMR/AMR)
- Mit BMI-Berechnung
- Interpretation des Körperfetts
- Elektroden aus gebürstetem Edelstahl
- 5 Aktivitätsgrade
- 10 Speicherplätze
- Überlastungsanzeige
- 180 kg Tragkraft
- Einteilung:
 - 0,1% für Körperfett, -wasser, Muskelanteil
 - 100 g für Gewicht, Knochenmasse
- Tap-On-Technik
- Abschaltautomatik
- Inkl. 2 x 3 V Batterie CR2032
- Produktabmessungen: 30 x 31,5 x 2,5 cm

5 Jahre Garantie

Versandereinheit: 4 Transportkarton: -
 EAN-Nr.: 42 11125 76430 2
 Art.-Nr.: 764.30

BG 21 Glass diagnostic scale

- Printing with "see-through"-effect
- Large LCD-display with multi display
- Digit size (weight display): 25 mm
- Body weight, body fat, body water, muscle ratio, bone mass, calory consumption (BMR/AMR)
- With BMI calculation
- Interpretation of body fat
- Electrodes of brushed stainless steel
- 5 activity levels
- 10 memories
- Overload indication
- 180 kg capacity
- Division:
 - 0,1% body fat, body water, muscle ratio
 - 100 g for body weight, bone mass
- Tap-on-technology
- Automatic switch off
- 2 x 3 V batteries CR2032 included
- Product measurements: 30 x 31.5 x 2.5 cm

5 year guarantee

Sales unit: 4 Shipping carton: -
 EAN no.: 42 11125 76430 2
 Item number: 764.30